

COMPETENT CELLS (H.M.M.: based on Hanahan method)

- STREAK a fresh plate the day before inoculation from parent stock (NOT previous batch of competent cells: glycerol stock only) with appropriate antibiotic using sterile technique and grow O/N 37°C.
 - XL1 or 2 Blue= TetR
 - TOP-10=StrepR (I plate without antibiotic for TOP-10)
 - Other strains: check with parent company (website). Most of our strains in the stock box are from Invitrogen or Stratagene
- INOCULATE a single colony into 5 ml **TYM** + antibiotic in a 14 ml Falcon 2059 tube; grow in shaker incubator at 37°C, 225 RPM.
- INOCULATE (in a.m.) 1 ml O/N culture into 500 ml **TYM** prewarmed to 37°C in a 2 l Erlenmeyer flask (NO antibiotic!!); grow in shaker incubator at 37°C, 225 RPM.
- GROW until OD550= 0.5 (approximately 2.5 hours); take sample to measure OD but keep the culture shaking and at 37°C
- COOL culture rapidly in ice water, 2' with swirling
- TRANSFER culture to pre-chilled 50 ml conical tubes
- SPIN 3K, 10', 4°C: NO BRAKE! (swinging bucket rotor –pre cool)
- REMOVE SUP, CAREFULLY tap pellet and GENTLY resuspend in 20 ml **Tfb1**/tube
 - try to get rid of clumps, but always be gentle
- INCUBATE ON ICE in cold room for 1.5 hours
- SPIN 3K, 10', 4°C: NO BRAKE!
- REMOVE SUP, CAREFULLY tap pellet and GENTLY resuspend in 1.5 ml **Tfb2**/tube
 - complete resuspension critical, be gentle
- ALIQUOT 200 µl to pre-labeled and pre-chilled eppendorfs
- FREEZE ON DRY ICE and STORE immediately at –80°C
- Calculate cfu/µg DNA using 3 serial 1/10 dilutions of 1 ng/µl pUC18

TYM: autoclave

2% Bacto-Tryptone
0.5% Bacto-Yeast
0.1 M NaCl
0.01 M MgCl₂

Tfb1: filter sterilize, store at 4°C

Make from autoclaved stocks

30 mM KOAc
50 mM MnCl₂ 4H₂O
100 mM KCl
10 mM CaCl₂
15% **w/v** glycerol
mix well, pH to 5.8 w/0.2 M acetic acid

Tfb2: filter sterilize, store at 4°C

Make from autoclaved stocks

10 mM Na-MOPS pH7.0

10 mM KCl

75 mM CaCl₂

15% w/v glycerol